

CIVCO's Needle Guidance Systems utilize a two-part system consisting of a reusable biopsy bracket and a disposable snap-on needle guide. The disposable needle guides snap onto the biopsy bracket, providing a fixed path for visualizing ultrasound-guided instrument placement procedures. Applications include tissue biopsy, fluid aspirations and catheter placement. Single-angle brackets offer imaging professionals a fixed angle needle path during puncture procedures. Brackets should be cleaned and disinfected according to your CIVCO User's Guide.

Ultra-Pro II™ Disposable Needle Guides

For use with VA57R-0375U, LO30E-075E, VA15R-050U, VA40R-035U, VA80R-035U Transducers

NEW ENHANCEMENTS • The Ultra-Pro II needle guide now incorporates more flexibility and advanced features to increase productivity during ultrasound-guided procedures. The needle guide directs instruments according to on-screen software guidelines. The Ultra-Pro II features a larger tab for improved quick-release function, allowing easy detachment of the needle from the transducer. Easy-to-read gauge sizes on the inserts makes it simple to identify and alter gauge sizes in a darkened ultrasound suite. The guide now features a larger funnel for instrument insertion and will accept a full range of sizes, including: 8.5FR, 14-23GA (19GA not available).



Ref. #	Product Description	Qty.
677-056	Non-sterile VA57R-0375U single-angle bracket (7cm)	24
610-699	Sterile Ultra-Pro II	24
610-579	Sterile Ultra-Pro II, 14 x 61cm (5.5" x 24") CIV-Flex cover (3D)	24
610-608	Sterile Ultra-Pro II, 14 x 91.5cm (5.5" x 36") telescopically-folded CIV-Flex cover (3D)	24
610-904	Sterile Ultra-Pro II, 14 x 147cm (5.5" x 58") telescopically-folded CIV-Flex cover (3D)	24

Multi-Pro 2000™ Disposable Needle Guides

For use with L038E-075E Transducer

The Multi-Pro 2000 provides a safe, efficient means of performing ultrasound-guided biopsy and instrument placement procedures. The Multi-Pro 2000 is widely used for ultrasound-guided fine-needle and core biopsy procedures, fluid aspirations and catheterization procedures. Inserts include 14, 16, 18, 20, 22 and 25 gauge sizes. Select kits include gel packet, transducer cover and colored elastic bands.



Ref. #	Product Description	Qty.
610-131	Sterile Multi-Pro 2000	24
610-133	Sterile Multi-Pro 2000, 14 x 61cm (5.5" x 24") CIV-Flex cover (3D)	24

Bi-Plane Transrectal Disposable Needle Guide

For use with UB10R-065U Transducer

This disposable needle guide's unique sliding design provides functionality for both biopsy and endocavity bi-plane scanning. The easy-to-use guide simply slides forward for use in biopsy procedures and retracts for transverse bi-plane scanning. Convenient single-use procedure kits include guide, covers, gel and bands. This endocavity guide accepts instruments from 16-18GA sizes.



Ref. #	Product Description	Qty.
677-057	Sterile transrectal needle guide	20
677-044	Sterile transrectal needle guide, 3.5 x 20cm (1.4" x 7.9") latex covers (2), 1.5 x 5.7cm (.6" x 2.3") latex fingercots	20

Reusable Transvaginal Needle Guide

For use with TV11R-055U Transducer

This stainless steel guide attaches to the probe by sliding into side locator grooves and can be secured by turning the thumbscrew into the closed position. To reuse the needle guide, follow standard cleaning, soaking and sterilization procedures as specified in your system's manual. The guide is sold as a non-sterile product. Sterile cover accessories are sold separately. This endocavity guide accepts instruments from 16-18GA sizes.



Ref. #	Product Description	Qty.
677-053	Non-sterile transvaginal needle guide, cleaning brush	1
610-667	Sterile 2 x 20cm (.8" x 7.9") latex cover, 3.5 x 20cm (1.4" x 7.9") latex cover, gel, bands and cleaning brushes (2)	20
610-324	Non-sterile endocavity cleaning brush, 17.8cm (7") 18GA	5

Transducers	CIV-Flex™ Latex-Free Loose-Fitting	NeoGuard™ Latex-Free Form-Fitting	Latex Form-Fitting	Polyethylene Latex-Free Loose-Fitting	Starter Kits Bracket Included	Replacement Kits Bracket Not Included
B096-035D	610-004		610-046	610-020		
I1038-050D			610-210	610-023, 610-217, 610-367		
I1038-050U			610-210	610-023, 610-217, 610-367		
I1072-050			610-210	610-023, 610-217, 610-367		
I1072-050D			610-210	610-023, 610-217, 610-367		
IT038-050D			610-210	610-023, 610-217, 610-367		
IT038-050U			610-210	610-023, 610-217, 610-367		
IT072-050			610-210	610-023, 610-217, 610-367		
L030E-075E	610-002, 610-542		610-044	610-020		610-579, 610-608, 610-699, 610-841
L038E-075E	610-002, 610-542		610-044	610-020		610-131, 610-133
L048-050D	610-004		610-046	610-020		
L072-035U	610-004		610-046	610-020		
L072-050D	610-002, 610-542		610-046	610-020		
L072-075D	610-002, 610-542		610-046	610-020		
LI10-035U	610-004		610-046	610-020		
LI23-035	610-004		610-046	610-020		
LI23-050	610-004		610-046	610-020		
S009-500	610-002, 610-542		610-044	610-020		
S011-375	610-001, 610-637		610-043	610-015		
S012-250	610-001, 610-637		610-043	610-015		
S012-250C	610-001, 610-637		610-043	610-015		
S111-375A	610-002, 610-542		610-044	610-020		
S112-250	610-002, 610-542		610-044	610-020		
S112-250CA	610-002, 610-542		610-044	610-020		
TV009-050U	610-006, 610-007	610-836, 610-843	610-010, 610-039, 610-075, 610-214			
TV11R-055U	610-006, 610-007	610-836, 610-843	610-010, 610-039, 610-075, 610-214		677-053	610-667
TV13R-050D	610-006, 610-007	610-836, 610-843	610-010, 610-039, 610-075, 610-214			
TV13R-050U	610-006, 610-007	610-836, 610-843	610-010, 610-039, 610-075, 610-214			
U072-050U	610-006, 610-007	610-836, 610-843	610-010, 610-039, 610-075, 610-214			
U072-050D	610-006, 610-007	610-836, 610-843	610-010, 610-039, 610-075, 610-214			
U096-050D	610-006, 610-007	610-836, 610-843	610-010, 610-039, 610-075, 610-214			
UB10R-065U	610-006, 610-007	610-836, 610-843	610-010, 610-039, 610-075, 610-214			677-044, 677-057
V15R-035D	610-002, 610-542		610-044	610-020		
V15R-035U	610-002, 610-542		610-044	610-020		
V15R-035U	610-002, 610-542		610-044	610-020		
V15R-050U	610-002, 610-542		610-044	610-020		610-579, 610-608, 610-699, 610-841
VA40R-035U	610-002, 610-542		610-044	610-020		610-579, 610-608, 610-699, 610-841
VA40R-050U	610-002, 610-542		610-044	610-020		
VAS7R-0375U	610-002, 610-542		610-044	610-020	677-056	610-579, 610-608, 610-699, 610-841
VA80R-035U	610-004		610-046	610-020		610-579, 610-608, 610-699, 610-841



Your Professional Connection for Quality Care

Worldwide Headquarters: 102 First Street South, Kalona, Iowa 52247-9589 USA
319/656-4447, 800/445-6741 • Fax: 319/656-4451, 877/329-2482 • www.civco.com

COPYRIGHT © 2005. ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL INSTRUMENTS CO. INC. CIV-FLEX, MULTI-PRO 2000, NEOGUARD, AND ULTRA-PRO II ARE ALL TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE THE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA 2003R-1623 REV. 3

to order call: 800.445.6741 • 319.656.4447