



ULTRASOUND IMAGING ACCESSORIES ECHOCARDIOGRAPHY



VOLUME 1
WWW.CIVCO.COM





ECHOCARDIOGRAPHY

Welcome to the 2010 CIVCO Medical Solutions Echocardiography catalog, part of our new Ultrasound Imaging Accessories catalog collection.

We are pleased you have chosen CIVCO's high-quality products and look forward to serving you. Our catalog will provide you with quick and easy access to the latest in echocardiography accessories for use with leading ultrasound equipment.

This year we divided our Ultrasound Imaging Accessories catalog into a collection of individualized specialty catalogs in order to better serve our customers. By offering products that meet clinicians' specific procedural needs, we hope to make the ordering process quicker and easier.

As a valued customer, CIVCO is committed to providing you with innovative and clinically beneficial products and services. Special care is taken for each and every product going out our door.

For over two decades, CIVCO has been providing imaging accessories to facilities around the world. Our company goal remains the same today, "Total Customer Satisfaction." If you are not completely satisfied with any CIVCO product, please contact one of our customer service specialists and they will be happy to assist you.

Thank you for your business. We look forward to serving you again in the near future.

*Visit WWW.CIVCO.COM for all of your
imaging accessory needs.*



Featured Solution



TEE Transducer Holder

For use with Acuson V5Ms, V7M, V7B, MPT7-4, Aloka UST-5293-5, UST-52116, UST-52104, UST-52110, UST-52111, GE Healthcare 6T, 6T-RS, 9T, 9T-RS, Philips microTEE, Omni III, S7-2t, S7-3t, T6207, X7-2t, SonoSite TEE Series and ZONARE T8-3 transducers

Visit WWW.CIVCO.COM/cardiology for more details or call 319.656.4447 for pricing

- Allows for stable and accurate positioning, while holding transducer in a fixed location for extended period of time
- Ergonomic design allows hands to be free, reducing the risk of repetitive stress injuries
- Attaches to most IV poles, tables and surgical rails
- Features a stable single-arm fixation with a single-point locking mechanism to quickly and easily lock the holder into position
- Reduces risk of stretching corners of patient's mouth over extended procedure time

REF. #	PRODUCT DESCRIPTION	QTY.
811-050	TEE Holder Kit with TEE Holder (1) Rail Adapter (1) and non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe and bite block (12) - Aloka, GE, Siemens, Sonosite, ZONARE	1
811-051	TEE Holder Kit with TEE Holder (1), Table Adapter (1), and non-sterile 4.6 tapered to 2.5 x 122 cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe, and bite block (12) - Aloka, GE, Siemens, Sonosite, ZONARE	1
811-052	TEE Holder Kit with TEE Holder(1) Rail Adapter(1) and non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel filled syringe and bite block(12) - Philips(Ref #8104)	1
811-053	TEE Holder Kit with TEE Holder(1), Table Adapter(1), and non-sterile 4.6 tapered to 2.5 x 122cm(1.8" tapered to 1"x48") CIV-Flex cover on applicator, gel filled syringe, and bite block(12) - Philips (Ref#8104)	1

Call 800.445.6741*
today to arrange
a FREE demo


* Inside USA Only



Prevention
Protection
Positioning

The TEE Transducer Holder and other CIVCO accessories are cost-effective solutions to these common problems:

- **Prevention** of cross-contamination for patient and staff safety.
- **Protection** of costly TEE transducer during transfer and storage.
- **Positioning** of TEE transducer for extended periods of time to improve staff ergonomics.



TEE Transducer Holder for the Operating Room

- Specifically for use in the OR and is compatible with most TEE probes
- Securely holds the handle of the probe
- Comes with a mounting clamp that attaches to a standard IV pole
- TEE Transducer Holder can be disinfected or sterilized by ethylene oxide, autoclave or MetriCide® and Cidex® cold chemical soaking

REF. #	PRODUCT DESCRIPTION	QTY.
610-729	TEE probe holder for the OR	1



Latex-Free Transesophageal Bite Blocks

- Tapered for tighter fit to TEE transducers
- Form-fitting cover provides defense against cross-contamination
- Bite blocks provide safe passage through patient's mouth
- Select kits include bite block and gel-filled syringe or gel packet

REF. #	PRODUCT DESCRIPTION	QTY.
610-160	Non-sterile E-Z Guard™ adult bite blocks without strap	24
610-749	Non-sterile E-Z Guard™ adult bite blocks with strap	24
610-979	Non-sterile latex-free E-Z Guard™ pediatric bite block with strap	24



ACUSON AcuNav™ Cover

For use with ACUSON AcuNav Ultrasound Catheter Family and GE Healthcare ICE Catheter Connector

- Provides latex-free barrier for the AcuNav catheter connector, handle and cable
- Accordion-fold allows for easy application
- Covers protect the integrity of the equipment & reduce risks of cross-contamination

REF. #	PRODUCT DESCRIPTION	QTY.
612-113	Sterile 15.2 tapered to 7.6 x 183cm (6" tapered to 3" x 72") polyethylene accordion-folded cover	12
612-227	Sterile 15.2 tapered to 7.6 x 183cm (6" tapered to 3" x 72") polyethylene accordion-folded cover	Demo Kit



Latex-Free CIV-Flex™ Adult TEE Cover Kits

- Tapered for tighter fit to TEE transducers
- Form-fitting cover provides defense against cross-contamination
- Protects TEE transducers
- Bite blocks provide safe passage through patient's mouth
- Select kits include bite block and gel-filled syringe or gel packet

REF. #	PRODUCT DESCRIPTION	QTY.
610-933	Sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel packet	12
610-840	Non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe, bite block	12
610-860	Non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator	12
610-1013	Sterile cover 3.8 tapered to 1.7 x 91.4cm (1.5" tapered to .7" x 36") CIV-Flex cover and gel packet	12
610-1014	Non-sterile cover 3.8 tapered to 1.7 x 91.4cm (1.5" tapered to .7" x 36") CIV-Flex cover, gel-filled syringe, and bite block	12



3M Red Dot™ ECG Electrodes

- Reliable traces, secure adhesion and comfortable fit for your patient
- Convenient, fast application and removal
- Latex-free and recommended for single-use application

REF. #	PRODUCT DESCRIPTION	QTY.
610-738	3M Red Dot Monitoring Electrode w/ clear-back, solid gel and snap (Ref # 2235)	50
610-737	3M Red Dot Monitoring Electrode w/ foam tape, solid gel and snap (Ref # 2259)	50
610-739	3M Red Dot Monitoring Electrode w/ soft cloth, solid gel and snap (Ref # 2238)	50



Lead-Lok® ECG Electrodes

- Pre-gelled for improved conductivity
- Offers foam backing and high-quality hypoallergenic adhesive
- Low impedance gel reduces the possibility of skin irritation
- Ensures reliable performance with clear, more stable ECG tracings
- Recommended for single-use application

REF. #	PRODUCT DESCRIPTION	QTY.
610-742	Lead-Lok ECG Monitoring Electrode w/ foam tape, gel and snap, 3 electrodes per strip (Ref # CIV 310 R)	150



Latex-Free Pediatric TEE Tip Guards

For use with Aloka UST-52111S and Philips S7-3t transducers

- Cylindrical plastisol protectors fit over the distal tip of a transesophageal transducer
- Protects transducer from damage during storage and transportation between patients
- Recommended for single-use

REF. #	PRODUCT DESCRIPTION	QTY.
610-1009	Non-sterile pediatric TEE tip guards	24



Visit WWW.CIVCO.COM for all of your imaging needs:

- Printing Supplies
- Electronic Storage
- Patient Leads & Cables
- TEE Leakage tester
- Transducer Covers
- Surgical Covers
- Gel & Gel Warmers
- Standoffs

CIVCO is a Complete Solutions Provider.



Latex-Free TEE Tip Guards

For use with Acuson/Siemens V5Ms, MPT7-4, Aloka UST-5293-5, UST-52104, UST-52116, Philips Omni III, S7-2 Omni, SonoSite TEE/8-3, Toshiba PEF-510MB and ZONARE T8-3 transducers

- Cylindrical plastisol protectors fit over the distal tip of a transesophageal transducer
- Protects transducer from damage during storage and transportation between patients
- Recommended for single-use
- Latex-Free tip guard

REF. #	PRODUCT DESCRIPTION	QTY.
610-945 (08265498)	Non-sterile TEE tip guard (yellow)	24



Philips X7-2 Latex-Free TEE Tip Guard

For use with Philips X7-2t transducer

- Cylindrical plastisol protectors fit over the distal tip of a transesophageal transducer
- Protects transducer from damage during storage and transportation between patients
- Recommended for single-use
- Latex-Free Tip Guard

REF. #	PRODUCT DESCRIPTION	QTY.
667-094	Non-sterile TEE tip guard-Philips X7-2t (gray) (Ref. # 453561326791)	24



NEW!

Latex-Free TEE Tip Guards

For use with Philips microTEE and S8-3t transducers

- Single-use cylindrical foam protectors fit over the distal tip of pediatric & adult transesophageal transducers
- Protects transducer from damage during storage and transportation between patients
- Recommended for single-use
- Latex-Free Tip Guard

REF. #	PRODUCT DESCRIPTION	QTY.
610-1110	Non-sterile TEE tip guard-Philips microTEE (white) (Ref. # 453561429541)	24
610-1111	Non-sterile TEE tip guard-Philips S8-3t (white) (Ref. # 453561409462)	24



Latex-Free Foam Tip Guards

- Single-use, cylindrical foam protectors
- Fits over distal tip of laparoscopic & pediatric TEE transducers
- Protects from damage during storage and transportation between patients
- Recommended for single-use
- Latex-Free Tip Guard

REF. #	PRODUCT DESCRIPTION	QTY.
610-748	Non-sterile foam tip guard	24

Vist WWW.CIVCO.COM/Cardiology or call a Customer Service Specialist at 319.656.4447 for additional product details.



TEE Fume-Free Soaking System

- Protects against tube and probe damage
- Packaged with two tubs, one for high-level disinfection and the other for initial rinse
- Meets all OSHA and JCAHO guidelines
- Dimensions: 48" x 14" x 15"
- Contact CIVCO at 319.656.4447 to order replacement filters

REF. #	PRODUCT DESCRIPTION	QTY.
610-850	Soaking station for TEE transducer, 110/220V (Ref. # G14TC-3)	1
610-851	Wall mounted storage tube (Ref. # AC-SUTC)	1
610-1100	GUS vapor control system filter, 110/60HZ (Ref.# F141)	1



TEE Probe Disinfection & Storage System

- Holds two transesophageal probes
- Securely houses and protects probes from physical damage
- Convenient for soaking processes
- Compatible with Cidex, MetriCide, and Omnicide & holds up to 1.4 liters
- Wall support dimensions: 67 x 28 x 11cm (26.5" x 11" x 4.25")

REF. #	PRODUCT DESCRIPTION	QTY.
610-724	TEE wall support	1
610-725	Storage container	1
610-723	Disinfection container	1



TEE Probe Storage Rack

- Provides convenient and safe storage
- Securely house and protects probes from physical damage
- Curved stain-relief prevents transducer cable from damage
- Supports a single probe
- Wall mounting hardware is not included

REF. #	PRODUCT DESCRIPTION	QTY.
610-1046	TEE transducer storage rack (Ref# AC-TCRK)	1



Transducer Storage Racks & Stabilizers

- Provides a safe, convenient storage location for endocavity and general purpose transducers
- Sturdy, metal rack supports up to four transducers and four stabilizers
- Wall mounting hardware included
- Transducer stabilizer sold separately

REF. #	PRODUCT DESCRIPTION	QTY.
610-942	Transducer storage rack - for use with Philips HD11, HDI 4000, HDI 5000, iU22 and iE33, and Siemens Antares transducers	1
644-070	Transducer storage rack - for use with Hitachi transducers	1
742-367	Transducer storage rack - for use with GE Healthcare and Siemens transducers	1
610-957	Transducer stabilizer - fits all transducer storage racks	1

COMPLETE SOLUTIONS PROVIDER

For more than 25 years, CIVCO has focused on designing, manufacturing and distributing products that improve patient care. Today our products help enhance life for patients undergoing ultrasound, radiation oncology, interventional radiology and urology services.

Contact us today to learn how you and your patients can benefit from a CIVCO solution.



For international distribution, visit CIVCO Worldwide at WWW.CIVCO.COM

Worldwide Headquarters: 102 First Street South, Kalona, Iowa 52247-9589 USA
Phone: 319.656.4447 | 800.445.6741 | Fax: 319.656.4451 | 877.329.2482
WWW.CIVCO.COM

COPYRIGHT © 2010. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIV-FLEX IS THE TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE THE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. PRINTED IN USA. 2010R-2479