



*Needle Guidance*



*Transducer Storage*



*In-Plane Needle Guidance*



*Equipment Protection*



# SIEMENS

Ultrasound Imaging Accessories • Volume 6

[WWW.CIVCO.COM](http://WWW.CIVCO.COM)

CIVCO's Needle Guidance Systems utilize a TWO-PART SYSTEM consisting of a REUSABLE bracket and a DISPOSABLE snap-on needle guide. The disposable needle guides snap onto the biopsy bracket, providing a fixed path for visualizing instrument placement on ultrasound-guided placement procedures. Each starter kit contains 5 sterile procedure kits including: needle guide, telescopically-folded CIV-Flex™ cover, gel packet and colored elastic band. Replacement kits are available from CIVCO. Brackets should be cleaned and disinfected according to your CIVCO User's Guide.

Each starter kit contains 5 sterile procedure kits including: needle guide, telescopically-folded CIV-Flex™ transducer cover, cleaning brushes, gel packet and colored elastic bands. U.S. Patent No. 5,941,889.

## Ultra-Pro II™ Reusable Starter Kits

### SG-1 Series: C5-2 • CX5-2

5.4, 9cm Depth



Ultra-Pro II™ style bracket  
**04838723**  
Replacement Kit 610-608  
Lock Pin Kit 610-1008

### SG-2 Series: PH4-1 • 4V1 (S2000)

5, 10cm Depth



Ultra-Pro II™ style bracket  
**07472330 (07471134)**  
Replacement Kit 610-608  
Lock Pin Kit 610-1007

### SG-3 Series: 9L4 VF7-3 • VF10-5 • VFX9-4

2.4, 4cm Depth



Ultra-Pro II™ style bracket  
**04838749**  
Replacement Kit 610-608  
Lock Pin Kit 610-1008

### SG-4 Series: 14L5 VF13-5 • VFX13-5

2.2cm Depth



Ultra-Pro II™ style bracket  
**04838756**  
Replacement Kit 610-608

### SG-5 Series: CH6-2

10cm Depth



Ultra-Pro II™ style bracket  
**07472322 (07471159)**  
Replacement Kit 610-608

### CH4-1 • CH5-2

4, 8cm Depth



Ultra-Pro II™ style bracket  
**07472314**  
Replacement Kit 610-609  
Lock Pin Kit 610-1007

### C5F1 • C7F2 • 7CF2

5, 10cm Depth



Ultra-Pro II™ style bracket  
**07843159**  
Replacement Kit 610-609  
Lock Pin Kit 610-1008

### 5C2g • 6C2

10cm Depth



Ultra-Pro II™ style bracket  
**08256947**  
Replacement Kit 610-608

### 4V1

5, 10cm Depth



Ultra-Pro II™ style bracket  
**08256950 (08256952)**  
Replacement Kit 610-608  
Lock Pin Kit 610-1007

### 4C1 • 6C1 HD

4, 8cm Depth



Ultra-Pro II™ style bracket  
**08260150 (08260151)**  
Replacement Kit 610-609  
Lock Pin Kit 610-1007

### L5 • L7 • 6L3 • 8L5

2.5cm Depth



Ultra-Pro II™ style bracket  
**612-085**  
Replacement Kit 610-608

### C3

11.4cm Depth



Ultra-Pro II™ style bracket  
**612-098**  
Replacement Kit 610-608

### V4

8.5cm Depth



Ultra-Pro II™ style bracket  
**612-099**  
Replacement Kit 610-608

### 17L5 HD

2.1, 5.4cm Depth



Ultra-Pro II™ style bracket  
**10037928 (10038223)**  
Replacement Kit 610-608  
Lock Pin Kit 610-1007

### 18L6 HD

2.1, 5.4cm Depth



Ultra-Pro II™ style bracket  
**10044200 (10044759)**  
Replacement Kit 610-608  
Lock Pin Kit 610-1007

### Multi-Angle Lock Pin



Replacement Kit:  
**610-1007 • 610-1008**

## Ultra-Pro II™ Disposable Needle Guides



- Directs devices according to on-screen software guidelines
- Large quick-release tab for easy detachment of needle from transducer
- Large funnel for device insertion
- Easy-to-read gauge size inserts
- Accepts 8.5FR & 14-23GA instruments (19GA not available)
- Designed for use with Ultra-Pro™ or Ultra-Pro II™ transducer brackets
- U.S. Patent No. 5,941,889

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-699</b>	Sterile Ultra-Pro II guide	24
<b>610-579</b>	Sterile Ultra-Pro II with 14 x 61cm (5.5" x 24") CIV-Flex cover (3D)	24
<b>610-608</b> <b>(07293074)</b>	Sterile Ultra-Pro II with 14 x 91.5cm (5.5" x 36") telescopically-folded CIV-Flex cover (3D)	24
<b>610-609</b>	Sterile Ultra-Pro II with 17.8 x 91.5cm (7" x 36") telescopically-folded CIV-Flex cover (3D)	24
<b>610-904</b>	Sterile Ultra-Pro II with 14 x 147cm (5.5" x 58") telescopically-folded CIV-Flex cover (3D)	24
<b>610-917</b>	Sterile Ultra-Pro II with 17.8 x 147cm (7" x 58") telescopically-folded CIV-Flex cover (3D)	24

## Ultra-Pro 3™ Disposable Needle Guides



- Zero needle resistance for smooth movement
- Quick-release function
- Accepts 11-22GA instruments (19, 21GA not available)
- Designed for use with Ultra-Pro or Ultra-Pro II transducer brackets
- U.S. Patent Nos. 5,941,889 and 7,909,815

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-918</b>	Sterile Ultra-Pro 3 guide	24
<b>610-901</b>	Sterile Ultra-Pro 3 with 14 x 91.5cm (5.5" x 36") telescopically-folded CIV-Flex cover (3D)	24
<b>610-919</b>	Sterile Ultra-Pro 3 with 17.8 x 91.5cm (7" x 36") telescopically-folded CIV-Flex cover (3D)	24

## Universal Needle Guide™ (UNG) Reusable Starter Kits



UNG™ needle guide

To order call Siemens at 800.228.4128, option 2

TRANSDUCER	ULTRASOUND SYSTEM
<b>3.5C40H</b>	ACUSON ANTARES
<b>5.0C50+</b>	SONOLINE G50, SONOLINE G60S
<b>5.0L45</b>	SONOLINE G50, SONOLINE G60S, ACUSON X500
<b>7.5L40</b>	ACUSON ANTARES
<b>7.5L70</b>	SONOLINE G50, SONOLINE G60S
<b>C5-2</b>	SONOLINE G20, SONOLINE G50
<b>C6-2</b>	SONOLINE G60S, ACUSON X500
<b>L10-5</b>	SONOLINE G20, SONOLINE G50, SONOLINE G60S, ACUSON X500

## Quik-Clip™ Reusable Starter Kits

Each starter kit contains 5 sterile procedure kits including: needle guide, telescopically-folded CIV-Flex™ transducer cover, cleaning brushes, gel packet and bands.

6L3 • 8L • L7 • L5



Quik-Clip™ style bracket  
**08235936**  
Replacement Kits: 612-062

To order call Siemens at  
800.228.4128, option 2



### Quik-Clip™ Disposable Replacement Kits

For use with 4V2, 5C2, 6L3, 8L5, C3, C366, C5, C544, L5, L582, L7, V2, and V4 transducers

- Plastic inserts for use with stainless steel Quik-Clip needle guide brackets
- Accepts 8FR, and 14, 16, 18, 20, and 22GA instruments
- Procedure kits include inserts, transducer cover, gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>612-062</b>	Sterile Quik-Clip Guide (small)	20
<b>612-058</b>	Sterile Quik-Clip for 4V2, 5C2, 6L3, 8L5, C3, C366, C5, C544, L5, L582, L7, V2, V4 barrel inserts, 5 x 30cm (2" x 11.8") latex cover	20
<b>612-068</b>	Sterile Quik-Clip for 4V2, 5C2, 6L3, 8L5, C3, C366, C5, C544, L5, L582, L7, V2, V4 barrel inserts, 14 x 61cm (5.5" x 24") CIV-Flex cover (3D)	20

## Quik-Clip™ Reusable Starter Kits

Each starter kit contains 5 sterile procedure kits including: needle guide, telescopically-folded CIV-Flex™ transducer cover, cleaning brushes, gel packet and bands.

8C4 • C7



Quik-Clip™ style bracket  
**08251375**  
Replacement Kits: 612-089

To order call Siemens at  
800.228.4128, option 2



### Quik-Clip™ Disposable Barrel Needle Guides

For use with 8C4 and C7 transducers

- Quik-Clip needle guide barrel inserts attach securely to reusable transducer bracket
- Barrel inserts provide accurate instrument placement for biopsy procedures
- Accepts 8FR and 14, 16, 18, 20, and 22 and 25GA instruments

REF. #	PRODUCT DESCRIPTION	QTY.
<b>612-089</b>	Sterile Quik-Clip Guide (large)	20



### Multi-Pro 2000™ Disposable Replacement Kits

For use with 3.5PL28, 7.5PL13, P8-3, S228, V228, S328, and V328 transducers

- Safe and efficient means for performing ultrasound-guided biopsies
- Includes inserts for 14, 16, 18, 20, 22 and 25GA instruments
- Select kits conveniently packaged with gel packet and bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-131</b>	Sterile Multi-Pro 2000 guide	24
<b>610-133</b>	Sterile Multi-Pro 2000 with 14 x 61cm (5.5" x 24") CIV-Flex cover (3D)	24
<b>610-134</b>	Sterile Multi-Pro 2000 with 17.8 x 61cm (7" x 24") CIV-Flex cover	24

Visit [WWW.CIVCO.COM](http://WWW.CIVCO.COM) for all of your imaging needs. CIVCO is your Solutions Provider.



**NEW!**

### Latex-Free Telescopically-Folded Super-Extended Length CIV-Flex™ Covers

- The new 96" cover provides complete transducer sterility from transducer array to connector
- Ultrasound covers improve patient infection control rates by preventing cross-contamination at minimal cost per procedure
- CIV-Flex material is a soft, pliable, latex-free material
- Sterile procedure kits include: transducer cover, sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-1119</b>	Sterile 10.2 x 244 cm (4" x 96") CIV-Flex cover	24
<b>610-1120</b>	Sterile 14 x 244 cm (5.5" x 96") CIV-Flex cover	24
<b>610-1121</b>	Sterile 17.8 x 244 cm (7" x 96") CIV-Flex cover	24



### Latex-Free Telescopically-Folded Extended-Length CIV-Flex Covers

- Extended lengths of 36" and 58" offer additional sterile protection
- Select covers offer 3D box end
- Sterile procedure kits include: transducer cover, sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-637</b>	Sterile 8.9 x 91.5cm (3.5" x 36") CIV-Flex cover	24
<b>610-1000</b>	Sterile 10.2 x 147cm (4" x 58") CIV-Flex cover	24
<b>610-1064</b>	Sterile 10.2 x 147cm (4" x 58") CIV-Flex cover (3D) (no gel)	24
<b>610-542</b>	Sterile 14 x 91.5cm (5.5" x 36") CIV-Flex cover (3D)	24
<b>610-575</b>	Sterile 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)	24
<b>610-1010</b>	Sterile 17.8 x 147cm (7" x 58") CIV-Flex cover (3D)	24



### Latex-Free Flat-Folded CIV-Flex Covers

- Offers sterile protection when performing routine puncture and drainage procedures
- Sterile procedure kits include: transducer cover, sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-323</b>	Sterile 6 x 91.5cm (2.5" x 36") CIV-Flex cover	24
<b>610-001</b>	Sterile 8.9 x 61cm (3.5" x 24") CIV-Flex cover	24
<b>610-002</b>	Sterile 14 x 61cm (5.5" x 24") CIV-Flex cover (3D)	24
<b>610-004 (04510947)</b>	Sterile 17.8 x 61cm (7" x 24") CIV-Flex cover	24
<b>610-005</b>	Sterile 22.9 x 61cm (9" x 24") CIV-Flex cover	24



### Latex-Free Bulk Packaged CIV-Flex Covers

- Bulk packaged in easy-to-use box
- Offers added protection against cross-contamination between routine scanning procedures
- Non-sterile colored elastic bands included

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-903</b>	Non-sterile 14 x 30cm (5.5" x 11.8") CIV-Flex cover (3D)	50
<b>610-362</b>	Non-sterile 14 x 61cm (5.5" x 24") CIV-Flex cover (3D)	50



### Latex-Free Polyethylene Covers

- Safe means for reducing cross-contamination
- Sterile procedure kits include: transducer cover, sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-013	Sterile 5 x 61cm (2" x 24") flat-folded polyethylene cover	24
610-015	Sterile 10.2 x 61cm (4" x 24") flat-folded polyethylene cover (3D)	24
610-1050	Sterile 15.2 x 25.4cm (6" x 10") flat-folded polyethylene cover (no gel)	24
610-208	Sterile 15.2 x 83.8cm (6" x 33") flat-folded polyethylene cover	24
610-989	Sterile 15.2 x 147cm (6" x 58") telescopically-folded polyethylene cover	24
610-020	Sterile 17.8 x 61cm (7" x 24") polyethylene cover	24
610-022	Sterile 22.9 x 61cm (9" x 24") polyethylene cover	24



### Latex Covers

- Rolled for easy application to a variety of transducers
- Offers sterile coverage during general puncture and drainage procedures
- Sterile procedure kits include: transducer cover, sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-042	Sterile 4 x 30cm (1.6" x 11.8") latex cover	24
610-043	Sterile 5 x 30cm (2" x 11.8") latex cover	24
610-044	Sterile 8 x 45cm (3.2" x 18") latex cover	24
610-046	Sterile 10 x 45cm (3.9" x 18") latex cover	24



### Transducer Storage Rack & Stabilizer

*For use with Siemens Antares transducers*

- Provides a safe, convenient storage location for endocavity and general purpose transducers
- Sturdy, metal rack supports up to four transducers and four stabilizers
- Wall mounting hardware included
- Transducer stabilizer sold separately

REF. #	PRODUCT DESCRIPTION	QTY.
610-942	Transducer storage rack	1
610-957	Transducer stabilizer	1












### Transducer Storage Rack & Stabilizer

*For use with Siemens transducers*

- Provides a safe, convenient storage location for endocavity and general purpose transducers
- Sturdy, metal rack supports up to four transducers and four stabilizers
- Wall mounting hardware included
- Transducer stabilizer sold separately

REF. #	PRODUCT DESCRIPTION	QTY.
742-367	Transducer storage rack	1
610-957	Transducer stabilizer	1

## Disposable Endocavity Needle Guides









	Transducer Compatibility	Product Description	REF. #	Qty.	Gauge Size	Guide Channel Length
	EV-8C4	Sterile, Disposable	676-154 10136306	24	16-18	16cm
	EC9-4 (for use on ACUSON Antares™ and S2000™ ultrasound systems)	Sterile, Disposable	610-842 610-711	24	16-18	16cm
	EC9-4 and EV9-4 transducers on SONOLINE G40™, ACUSON X150™ and X300™ ultrasound systems	Sterile, Disposable	676-132 10041822	24	16-18	14cm
	EC7, EC10C5, EC9-4 and EV9-4 transducers (EC9-4 on SONOLINE G20™, G50™, G60 S™; EC9-4 and EV9-4 on ACUSON X500™ ultrasound systems)	Sterile, Disposable	612-061 612-051	24	16-18	12cm
	6.5EC10	Sterile, Disposable	657-014 610-543	24	16-18	12cm
	BP9-4	Sterile, Disposable	676-150 676-144	24	16-18	11cm
	BE9-4 (short)	Sterile, Disposable	676-112 676-110	24	16-18	12cm
	BE9-4 (long)	Sterile, Disposable	610-590 610-693 610-588	24	16-18	14cm
	V714S and I714S	Sterile, Disposable	612-045 612-035	24	16-18	12cm



*CIVCO Medical Solutions is the leader in infection prevention and offers numerous single-use products to protect patients and help hospitals reduce hospital-acquired infections resulting from cross-contamination.*

**Keep it Sterile with CIVCO.**

## Reusable Endocavity Needle Guides

	Transducer Compatibility	Product Description	REF. #	Qty.	Gauge Size	Guide Channel Length
	EV-8C4 and EVC8	Non-sterile, Reusable	08245962	1	16-18	14cm
	EC9-4 and EV9-4 transducers on the SONOLINE G40™, ACUSON X150™ and X300™ ultrasound systems	Non-sterile, Reusable	10041823	1	16-18	14cm
	EC7, EC10C5, 6.5EC10, EC9-4 and EV9-4 transducers (EC9-4 on SONOLINE G20™, G50™, G60 S™; EC9-4 and EV9-4 on ACUSON X500™ ultrasound systems)	Non-sterile, Reusable	610-666	1	16-18	12cm
	EV9F4 and 9EVF-4	Non-sterile, Reusable	10031556	1	16-22	18cm
	6.5EV13, 6.5EV13+, 6.5EV13S and EV9-4 transducers (EV9-4 on SONOLINE G20™, G50™, G60 S™ ultrasound systems)	Non-sterile, Reusable	04392143	1	16-22	19cm
	BP9-4 U.S. Patent Pending	Non-sterile, Reusable	10563772	1	18 and smaller	11.2cm
	BE9-4 (short)	Non-sterile, Reusable	08651387	1	16-18	12cm
	BE9-4 (long)	Non-sterile, Reusable	10035994	1	16-18	14cm

Note: All reusable needle guides ship with cleaning brushes (5).



### Cleaning Brushes

- Designed for cleaning reusable CIVCO needle guide cannulas
- Recommended for single-use with any reusable endocavity needle guide
- Manufactured from stainless steel wire
- Brushes offer 1/2" finger ring

REF. #	PRODUCT DESCRIPTION	QTY.
610-324	Endocavity needle guide cleaning brush, 22.9cm (9"), 18GA	5
610-951	Endocavity needle guide cleaning brush, 22.9cm (9"), 18GA	200



### Endocavity Thumbscrew Replacement Kit

For use with Siemens reusable needle guides

- New stainless steel thumbscrew closure replacement kits now available
- Designed and recommended for use with probe specific needle guides
- Contact CIVCO for compatibility details

REF. #	PRODUCT DESCRIPTION	QTY.
610-1025	Non-sterile endocavity thumbscrew (Ref. #079-400)	5



### Latex-Free NeoGuard™ Covers

- Rolled, latex-free material designed for ultrasound use
- Recommended alternative for latex-sensitive patients and professionals
- Covers include colored elastic bands
- Sterile gel packets available separately

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-843</b>	Sterile 2.6 x 30cm (1" x 11.8") NeoGuard cover	24
<b>610-1038</b>	Sterile 2.6 x 30cm (1" x 11.8") NeoGuard cover	50
<b>610-836</b>	Non-sterile 2.6 x 30cm (1" x 11.8") NeoGuard cover	24
<b>610-1039</b>	Non-sterile 2.6 x 30cm (1" x 11.8") NeoGuard cover	50
<b>676-117 (10031547)</b>	Non-sterile 2.6 x 30.5cm (1" x 12") NeoGuard cover for use with Endo-Pil onG20 ultrasound system	24
<b>610-844</b>	Sterile 4 x 30cm (1.6" x 11.8") NeoGuard cover	24
<b>610-838</b>	Non-sterile 4 x 30cm (1.6" x 11.8") NeoGuard cover	24
<b>610-1126</b>	<b>NEW</b> • Sterile 2 x 30cm (.8" x 11.8") NeoGuard cover	24
<b>610-1127</b>	<b>NEW</b> • Non-sterile 2 x 30cm (.8" x 11.8") NeoGuard cover	24



### Latex-Free Eclipse® Cover

- Now available for both 2D and 3D/4D ultrasound transducers
- Latex-free, pliable, polyethylene material
- Pre-gelled with Aquasonic® 100 transmission gel
- Reduces concerns of cross-contamination and protects the transducer during procedures

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-792</b>	Non-sterile 6.3 tapered to 4.6 x 24.1cm (2.5" tapered to 1.8" x 9.5") polyethylene cover with gel	100
<b>610-1102</b>	<b>NEW</b> • Non-sterile 8.3 tapered to 4.1 x 24.1cm (3.3" tapered to 1.6" x 9.5") polyethylene cover with gel (Ref. # 38-03)	100



### Latex-Free CIV-Flex™ Covers

- Soft, pliable material provides comfort
- Acoustically transparent, offering distortion-free imaging
- Slides easily over transducers
- Recommended alternative for latex-sensitive patients and professionals
- Sterile cover kits include gel and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-006</b>	Sterile 11.9 tapered to 3.8 x 61cm (4.7" tapered to 1.5" x 24") flat-folded CIV-Flex cover, gel	24
<b>610-007</b>	Non-sterile 11.9 tapered to 3.8 x 61cm (4.7" tapered to 1.5" x 24") flat-folded CIV-Flex cover (no gel)	24

**Visit [WWW.CIVCO.COM/SIEMENS](http://WWW.CIVCO.COM/SIEMENS) to find additional products recommended for use with your SIEMENS Ultrasound System.**



### Latex PROcovers™

- High-quality, custom-sized
- Provides defense against cross-contamination
- Includes colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-075</b>	Sterile 2 x 20cm (.8" x 7.9") latex cover	24
<b>610-039</b>	Non-sterile 2 x 20cm (.8" x 7.9") latex cover	24
<b>610-212</b>	Sterile 2 x 30cm (.8" x 11.8") latex cover	24
<b>610-076</b>	Sterile 2.6 x 20cm (1" x 7.9") latex cover	24
<b>610-040</b>	Non-sterile 2.6 x 20cm (1" x 7.9") latex cover	24
<b>610-213</b>	Sterile 2.6 x 30cm (1" x 11.8") latex cover	24
<b>610-101</b>	Non-sterile 2.6 x 30cm (1" x 11.8") latex cover	24
<b>610-042</b>	Sterile 4 x 30cm (1.6" x 11.8") latex cover	24
<b>610-041</b>	Non-sterile 4 x 30cm (1.6" x 11.8") latex cover	24



### Tapered Latex Covers

- Protects against cross- contamination
- Provides form-fitting tip
- Offers easy application to larger handled transducers
- Sterile procedure kit includes transducer cover, gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-985</b>	Sterile 3.2 tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") latex cover	24
<b>08990426</b>	Non-sterile 3.2 tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") latex cover with polyethylene cord cover, gel and elastic bands	10
<b>610-795</b>	Non-sterile 3.2 tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") latex cover	24
<b>08990418</b>	Non-sterile 3.2 tapered to 1.9 x 32cm (1.3" tapered to .75" x 12.6") latex cover	100



### 3.5 x 20cm Bulk Packaged Latex Covers

- As low as 19 cents per cover
- Packaged in protective light-sensitive bags for extended shelf life
- Provides protection against cross-contamination
- Colored elastic bands sold separately

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-249</b>	Non-sterile 3.5 x 20cm (1.4" x 7.9") latex cover	100
<b>610-262</b>	Non-sterile 3.5 x 20cm (1.4" x 7.9") latex cover	500



### 3.5 x 20cm Latex Covers

- Conveniently packaged in individual perforated packets
- Looser fitting rolled latex cover
- Provides protection from cross-contamination
- Elastic bands included with each cover

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-214</b>	Sterile 3.5 x 20cm (1.4" x 7.9") latex cover	50
<b>610-011</b>	Non-sterile 3.5 x 20cm (1.4" x 7.9") latex cover	24
<b>610-010</b>	Non-sterile 3.5 x 20cm (1.4" x 7.9") latex cover	50



## Latex-Free Endocavity Balloon

For use with Siemens Endo-P II transducer

- Provides small insertion profile to minimize patient discomfort
- Features adjustable inflation control to maximize user control over prostate shape and height
- Fits snugly and enhances ultrasound image
- U.S. Patent No. 6,126,607

REF. #	PRODUCT DESCRIPTION	QTY.
610-898	Non-sterile 2 x 14cm (.8" x 5.5") Endocavity Balloon	10



## Latex Waterpath Balloon

For use with Siemens Endo-P II transducer

- Designed with an integrated tube to allow easy filling and deventilation of balloon
- Designed for enhanced acoustic coupling and waterpath standoff to rectal wall
- Serves as a hygienic transducer cover
- Product contains latex

REF. #	PRODUCT DESCRIPTION	QTY.
610-827 (09704461)	Non-sterile 2 x 19cm (.8" x 7.5") latex balloon	50



## Endocavity Transducer Soaking Cup

- Accommodates two transducers, including larger 3D probes
- Recommended for use with glutaraldehyde cold chemical disinfecting/sterilizing solutions
- Holds up to 120 fl. oz. of disinfecting liquid
- Kit includes 12" x 5" polyethylene cup, lid, mounting hardware and contoured funnels for smaller transducers

REF. #	PRODUCT DESCRIPTION	QTY.
610-584	Starter kit with bracket, polyethylene cup, lid, smaller-white funnel, larger-blue funnel and mounting hardware	1
610-585	Replacement kit with cup, lid, smaller-white funnel (3) and larger-blue funnel (3)	3
610-586	Soaking cup stabilizer, rubber feet (5), assembly hardware	1

Use CIVCO's Transducer Search to find compatible needle guidance solutions, transducer covers and other imaging accessories available for use with your Siemens transducer and ultrasound system.

[WWW.CIVCO.COM/TRANSDUCERSEARCH](http://WWW.CIVCO.COM/TRANSDUCERSEARCH)





### Latex-Free Surgi™ Intraoperative Covers

- Polyethylene cover accordion-folded for easy application to transducers
- Extends up transducer and cord 2.4 meters (eight feet) for additional sterile protection
- Flat-folded cover also available
- Select covers include gel packets and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-023</b>	Sterile 15.2 tapered to 7.6 x 244cm (6" tapered to 3" x 96") accordion-folded cover (no gel)	12
<b>610-367</b>	Sterile 15.2 tapered to 7.6 x 244cm (6" tapered to 3" x 96") accordion-folded cover	12
<b>610-217</b>	Sterile 15.2 tapered to 7.6 x 244cm (6" tapered to 3" x 96") flat-folded cover	24



### Latex-Free Surgi-Tip™ Intraoperative Covers

- Provides snug fit to small footprint transducers
- Accordion or telescopically-folded for easy application
- Extended sterile coverage for transducer and cord
- Covers include: sterile gel packet, probe clips and colored elastic bands


REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-1072</b>	Sterile 15.2 x 244cm (6" x 96") accordion-folded cover with attached 3cm (1.2") latex-free NeoGuard tip with clips	12
<b>610-833</b>	Sterile 15.2 x 244cm (6" x 96") poly telescopically-folded cover with attached 3cm (1.2") latex-free NeoGuard tip with clips	12
<b>610-834</b>	Sterile 20.3 x 244cm (8" x 96") poly telescopically-folded cover with attached 6cm (2.4") latex-free NeoGuard tip with clips	12



### Intraoperative Non-Pyrogenic Transducer Covers

- Provides sterile, non-pyrogenic barrier
- Protects against cross-contamination
- Prevents transducer contact with patients suffering from Creutzfeldt-Jakob Disease (CJD)
- Latex-free, CIV-Flex™ cover includes sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-956</b>	Sterile non-pyrogenic 7.6 tapered to 4.1 x 147cm (3" tapered to 1.6 x 58") telescopically-folded CIV-Flex cover	24
<b>610-1107</b>	Sterile non-pyrogenic 8.9 x 147cm (3.5"x 58") telescopically-folded CIV-Flex cover	24



Visit [WWW.CIVCO.COM](http://WWW.CIVCO.COM) for all of your imaging needs.  
**CIVCO is your Complete Solutions Provider.**



### ACUSON AcuNav™ Covers

For use with ACUSON AcuNav Ultrasound Catheter Family Connector

- Provides latex-free barrier for AcuNav catheter connector, handle and cable
- Accordion-folded covers apply easily
- Covers protect the integrity of equipment and reduce risks of cross-contamination
- Individually packaged with tape strip and clips to secure the cover

REF. #	PRODUCT DESCRIPTION	QTY.
<b>612-113</b> <b>(08264920)</b>	Sterile 15.2 tapered to 7.6 x 183cm (6" tapered to 3" x 72") accordion-folded polyethylene cover	12




### Latex-Free Surgi-Boot™ Intraoperative Covers

For use with Siemens VF13-5SP transducers

- Extends field of sterile protection against cross-contamination
- Latex-free, form-fitting tip provides snug fit to small footprint transducers
- Individually packaged with sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-797</b> <b>(07693646)</b>	Sterile 15.2 x 244cm (6" x 96") poly telescoped cover with attached form-fitting latex-free NeoGuard tip	12



### Latex-Free Laparoscopic Cover

- Form-fitting tipped transducer cover fits standard 10mm laparoscopic transducers (requires 12mm trocar)
- Extends up transducer and cord 2.4 meters (eight feet) for additional sterile protection
- Individually packaged with surgical lubricant and probe clips

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-941</b> <b>(08247399)</b>	Sterile 15.2 x 244cm (6" x 96") poly telescoped cover with attached form-fitting latex-free tip with clips	12



### Latex-Free Intraoperative System Drapes & Cord Covers

- Cord cover and system drape provides wide sterile field
- Probe clips and colored elastic bands included with all cord covers
- Adhesive tape included with each system drape

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-1037</b>	Sterile 50.8 x 71.1 cm (20" x 28") polyethylene system drape, tape strips (4)	24
<b>610-025</b> <b>(07693653)</b>	Sterile 76 x 76cm (30" x 30") polyethylene system drape	24
<b>610-024</b>	Sterile 7.6 x 244cm (3" x 96") polyethylene cord cover	24



## TEE Transducer Holder

For use with Acuson V5Ms, V7M, V7B, MPT7-4 transducers

Visit [WWW.CIVCO.COM/cardiology](http://WWW.CIVCO.COM/cardiology) for more details

- Allows for stable and accurate positioning, while holding transducer in a fixed location for extended period of time
- Ergonomic design allows hands to be free, reducing the risk of repetitive stress injuries
- Attaches to most IV poles, tables and surgical rails
- Features a stable single-arm fixation with a single-point locking mechanism to quickly and easily lock the holder into position
- Reduces risk of stretching corners of patient's mouth over extended procedure time

REF. #	PRODUCT DESCRIPTION	QTY.
<b>811-050</b> <b>(811-054)</b>	TEE Holder Kit with TEE Holder (1) Rail Adapter (1) and non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe and bite block (12)	1
<b>811-051</b> <b>(811-055)</b>	TEE Holder Kit with TEE Holder (1), Table Adapter (1), and non-sterile 4.6 tapered to 2.5 x 122 cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe, and bite block (12)	1



## Latex-Free CIV-Flex™ TEE Covers & Bite Blocks

- Custom tapered design provides tighter fit to adult, small adult and pediatric TEE transducers
- Cover provides a defense against cross-contamination while protecting the transducer
- Form-fitting cover is preloaded on applicator for ease of application to TEE transducer
- Bite blocks provide a safe passage through the patient's mouth.

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-979</b>	Non-sterile latex-free E-Z Guard™ pediatric bite block with strap	24
<b>610-160</b>	Non-sterile latex-free E-Z Guard™ adult bite blocks without strap	24
<b>610-749</b>	Non-sterile latex-free E-Z Guard™ adult bite blocks with strap	24
<b>610-933</b>	Sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator and gel packet	12
<b>610-840</b>	Non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe and bite block	12
<b>610-860</b>	Non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator	12
<b>610-1013</b>	Sterile cover 3.8 tapered to 1.7 x 91.4cm (1.5" tapered to .7" x 36") CIV-Flex cover and gel packet	12
<b>610-1014</b>	Non-sterile cover 3.8 tapered to 1.7 x 91.4cm (1.5" tapered to .7" x 36") CIV-Flex cover, gel-filled syringe, and bite block	12



## Latex-Free TEE Tip Guards

For use with MPT7-4 and V5Ms transducers

- Cylindrical plastisol protectors fit over the distal tip of a transesophageal transducer
- Protects transducer from damage during storage and transportation between patients
- Recommended for single-use
- Latex-free tip guard

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-945</b> <b>(08265498)</b>	Non-sterile TEE tip guard (yellow)	24



### TEE Fume-Free Soaking System

- Protects against tube and probe damage
- Packaged with two tubes, one for high-level disinfection and the other for initial rinse
- Meets all OSHA and JCAHO guidelines
- Dimensions: 51" x 17" x 16"
- Weight: 55 lbs.
- Contact CIVCO at 319.656.4447 to order replacement filters

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-850</b>	Soaking station for TEE transducer, 110/220V (Ref. # G14-TC3)	1
<b>610-851</b>	Wall mounted storage tube (Ref. # AC-SUTC)	1
<b>610-1100</b>	GUS vapor control system filter, 110/60HZ (Ref.# F141)	1



### TEE Probe Disinfection & Storage System

- Holds two transesophageal probes
- Securely houses and protects probes from physical damage
- Convenient for soaking processes
- Compatible with Cidex, MetriCide, and Omnicide & holds up to 1.4 liters
- Wall support dimensions: 67 x 28 x 11cm (26.5" x 11" x 4.25")

REF. #	PRODUCT DESCRIPTION	QTY.
<b>610-724</b>	TEE wall support	1
<b>610-725</b>	Storage container	1
<b>610-723</b>	Disinfection container	1

BANDS & CLIPS	REF. #	QTY.
Non-sterile latex-free bands, 1.6cm (.63"), large, orange	610-689	100
Non-sterile probe clips, small	610-125	100
Non-sterile probe clips, VNAV	610-1033	25
RESPIRATION MONITOR ACCESSORIES	REF. #	QTY.
Nasal probe (1), disposable nose clips (5)	4000-0716-02	1
Non-sterile disposable nose clips	8000-1373-01	25
PATIENT POSITIONING	REF. #	QTY.
<b>NEW</b> • Non-sterile Kanguru Web Universal Base Unit (Ref. #RB2555)	610-1047	1
<b>NEW</b> • Non-sterile Kanguru Web Disposable Webs (Ref. #DW1000)	610-1048	10
GENERAL PURPOSE STANDOFFS	REF. #	QTY.
Disposable Aquaflex® gel pads (Ref. #04-02)	610-113	6
Acoustic standoff, 1 x 10 x 15cm (.4" x 4" x 6")	610-459	1
Acoustic standoff, 2 x 10 x 15cm (.8" x 4" x 6")	610-460	1
Sterile Hydroscan® couplant sheets, 10 x 15cm (4" x 6")	610-346	5

TRANSDUCER POSITIONING	REF. #	QTY.
Non-sterile positioning arm kit with positioning arm 32.4" (82.3cm) (Ref. #810-200), rail adapter-90 degree, 1.1" x .38" (2.8cm x 1cm) (Ref. #810-106), multi-purpose holder (Ref. #810-202) and sterile 25.4 tapered to 5.1 x 121.9cm (10" tapered to 2" x 48") telescopically-folded CIV-Flex cover with Quick-Mount™ perforation (5) (Ref. #610-932)	811-001	1
Non-sterile positioning arm kit with positioning arm 32.4" (82.3cm) (Ref. #810-200), table adapter 1" (2.5cm) (Ref. #810-100), multi-purpose holder (Ref. #810-202) and sterile 25.4 tapered to 5.1 x 121.9cm (10" tapered to 2" x 48") telescopically-folded CIV-Flex cover with Quick-Mount perforation (5) (Ref. #610-932)	811-002	1
Non-sterile positioning arm kit with positioning arm-extended length, 40.4" (102.6cm) (Ref. #810-205), rail adapter- 90 degree, 1.1" x .38" (2.8cm x 1cm) (Ref. #810-106), multi-purpose holder (Ref. #810-202) and sterile 25.4 tapered to 5.1 x 121.9cm (10" tapered to 2" x 48") telescopically-folded CIV-Flex cover with Quick-Mount perforation (5) (Ref. #610-932)	811-003	1
Non-sterile positioning arm kit with positioning arm-extended length, 40.4" (102.6cm) (Ref. #810-205), table adapter 1" (2.5cm) (Ref. #810-100), multi-purpose holder (Ref. #810-202) and sterile 25.4 tapered to 5.1 x 121.9cm (10" tapered to 2" x 48") telescopically-folded CIV-Flex cover with Quick-Mount perforation (5) (Ref. #610-932)	811-004	1
Sterile 10" tapered to 2" x 48" (25.4 tapered to 5.1 x 121.9cm) telescopically-folded CIV-Flex cover with Quick-Mount perforation	610-932	48
ENDOCAVITY ACCESSORIES	REF. #	QTY.
Non-sterile 2 x 7cm (.8" x 2.8") latex fingertip	610-038	100
Non-sterile 1.5 x 7cm (.6" x 2.8") latex fingertip	610-316	100
Non-sterile latex-free o-rings, 1.3cm (.5"), black	612-052	40
Non-sterile latex-free o-rings, 1.6cm (.63"), black	610-751	100
Non-sterile latex-free o-rings, 1.9cm (.75"), white	610-145	100
Non-sterile latex-free o-rings, 2.3cm (.92"), white	610-752	25
Non-sterile latex bands, .8cm (.31"), small, green	610-690	500
Non-sterile latex bands, .8cm (.31"), small, green	610-691	10,000
ENDOCAVITY SOAKING STATION	REF. #	QTY.
GUS vapor control system filter, 110/60HZ (Ref. #F10)	610-1101	1
Soaking station for one endocavity transducer, 110V/60Hz (Ref. #G10VP)	610-921	1
Soaking station for two endocavity transducers, 110/60Hz (Ref. #G14KA)	610-720	1
SCANNING KITS	REF. #	QTY.
Non-sterile scanning kit with extension tube, syringe and stopcock	610-224	5
Sterile extension tube 53.3cm (21")	610-256	10
Non-sterile scanning kit with syringe (1), extension tube 53.3cm (21") (1), 3-way stopcock (1), 2 x 20cm (.8" x 7.9") latex covers (24)	612-056	24
Non-sterile 3.5 x 20cm (1.4" x 7.9") latex cover with syringe (1), extension tube (1), stopcock (1) and a bag of bands (1)	610-058	48
PATIENT LEADS , CABLES & ELECTRODES	REF. #	QTY.
Clear-backed electrode with snap (Ref. #2235)	610-738	50
Foam monitoring electrode with snap (Ref. #2259)	610-737	50
Soft cloth electrode with snap (Ref. #2238)	610-739	50
Resting electrode - adult (Ref. #CIV 310 R)	610-742	150
TEE POSITIONING AND STORAGE	REF. #	QTY.
TEE probe holder for the OR	610-729	1
TEE transducer storage rack (Ref. #AC-TCRK)	610-1046	1

ULTRASOUND LOTION	REF. #	QTY.
Non-sterile Polysonic Ultrasound Lotion, 250ml bottles (Ref. #21-08)	610-1139	12
Non-sterile Polysonic Ultrasound Lotion, 3.8 liters with dispenser (Ref. #21-28)	610-1140	1
ULTRASOUND GEL	REF. #	QTY.
Aquagel® 5 oz. flip-top tube (Ref. #57-05)	610-920	20
Aquasonic® 100, 20 gram sterile gel packets (Ref. #01-01)	610-117	48
Aquasonic® 100, .25ml bottles (Ref. #01-08)	610-130	12
Aquasonic 100, 100-20 gram gel packettes per dispenser box (Ref. #01-20)	610-1122	4
Aquasonic 100, 5 liters, 1 dispenser bottle (Ref. #01-50)	610-116	1
Aquasonic Clear, 5 liters econopack, 2 dispenser bottles (Ref. #03-54)	610-455	4
Scan™, 237ml bottle (Ref. #11-08)	610-892	12
Scan, 1 gallon pump, 2 dispenser bottles (Ref. #11-28S)	610-123	4
Ultra/Phonic™, 20ml sterile gel packets (Ref. #NDC36-1000-20)	610-057	25
GEL WARMERS	REF. #	QTY.
Gel and lotion warmer, three-bottle 110V (Ref. #201-786)	610-993	1
NEW • Thermasonic™, single-bottle gel warmer-120V (Ref. #82-01)	610-1020	1
NEW • Thermasonic, three-bottle gel warmer-120V (Ref. #82-03)	610-1022	1
SKIN MARKING PENS	REF. #	QTY.
Non-sterile SQUEEZE-MARK® skin marking pen with ruler - medium tip	610-746	100
DISINFECTANTS	REF. #	QTY.
Cidex®, 14 day use, 1.25 gallon bottle (Ref. #735309)	610-640	4
Cidex, 14 day use test strips, 60 strips per bottle	610-803	2
Cidex Plus®, 28 day use, 1 gallon bottle (Ref. #2785)	610-641	4
Cidex Plus, 28 day use test strips, 60 strips per bottle	610-804	2
Cidex OPA®, 14 day use, 1 gallon bottle (Ref. #20390)	610-704	4
Cidex OPA test strips, 60 strips per bottle	610-714	2
MetriCide® 28, 28 day use, 1 quart bottle (Ref. #2805)	610-328	16
MetriTest® strips, 60 strips per bottle	610-805	2
Enzol®, 1 gallon bottle (Ref. #2252)	610-643	4
NEW • Glute-Out neutralizer, 2 oz. bottle (Ref. #NG-1GAL)	610-1045	24
Glute-Out neutralizer, .5oz bottles (Ref. #NG-QT)	610-1099	24
Spill kit with Glute-Out Neutralizer NG-5GAL bottle (2), impervious gown (2), organic respirator (1), waterproof shoes (2 pair), nitrile gloves (2 pair), safety goggles (1), absorbent mats (4), and bouffant caps (2) (Ref. #AC-SP1)	610-1098	1
DISINFECTANTS SPRAYS	REF. #	QTY.
DISPATCH® liquid spray, 22 oz. (Ref. #68967)	610-705	8
Protex™ disinfectant spray, 12 oz. (Ref. #42-12)	610-1079	12
T-Spray™ transducer spray, 250ml pump spray (Ref. #NDC36-3200-25)	610-154	12
T-Spray II™ transducer spray, 250ml pump spray (Ref. #NDC36-2102-25)	610-703	12
DISINFECTANTS TOWELETTES	REF. #	QTY.
Sani-Cloth HB® Germicidal Disposable Wipe, 15.3 x 17.1cm (6" x 6.75"), 160 sheets per dispenser (Ref. #Q08472)	610-991	12
Sani-Cloth HB Germicidal Disposable Wipe, 20.3 x 35.6cm (8" x 14"), 65 sheets per dispenser (Ref. #Q85484)	610-992	6
Sani-Cloth® PLUS Germicidal Dipsoable Wipe, (6" x 6.75") 160 sheets per dispenser (Ref. #Q89072)	610-943	12

Sani-Cloth PLUS Germicidal Disposable Wipe, (8" x 14") 65 sheets per dispenser (Ref. #Q85084)	610-944	6
<b>CLEANSING SPRAYS</b>	<b>REF. #</b>	<b>QTY.</b>
Transeptic® cleansing solution, 250ml pump spray (Ref. #09-25)	610-388	12
<b>CLEANSING TOWELETTES</b>	<b>REF. #</b>	<b>QTY.</b>
Towelettes, (7" x 6") 85 sheets per dispenser	610-128	6
<b>PRINTING SUPPLIES</b>	<b>REF. #</b>	<b>QTY.</b>
Ultrasound Print Frames	610-733	100
Chart insert for 1-4, 110mm thermal prints	610-735	200
Chart insert for 1-12 thermal prints	610-732	150
Top-binding print holder, 1-8 thermal prints (plastic)	610-734	100
<b>FIDUCIAL MARKERS</b>	<b>REF. #</b>	<b>QTY.</b>
Non-sterile PinPoint® multi-modality fiducial marker with general purpose sensor connection 1.3mm	610-1053	1
Non-sterile PinPoint multi-modality fiducial markers, 1.3mm diameter (Ref. #128)	610-1054	1
<b>OPTICAL DISKS &amp; CDs</b>	<b>REF. #</b>	<b>QTY.</b>
Sony 3.5" MO Media, 128MB MOD, 512 bytes/sector (Ref. #EDM-128C)	610-630	5
Sony 3.5" MO Media, 230MB MOD, 512 bytes/sector (Ref. #EDM-230C)	610-631	5
Sony 3.5" MO Media, 540MB MOD, 512 bytes/sector (Ref. #EDM-540C2)	610-632	5
Sony 3.5" MO Media, 640MB MOD, 2048 bytes/sector (Ref. #EDM-640C2)	610-633	5
Mediscrite 3.5" MO Media, 230MB MOD, 512 bytes/sector (Ref. #MO-R230MBX5)	610-894	5
Mediscrite 3.5" MO Media, 640MB MOD, 512 bytes/sector (Ref. #MO-R640MBX5)	610-895	5
<b>MITSUBISHI PRINT MEDIA</b>	<b>REF. #</b>	<b>QTY.</b>
KP-61S b/w super grade paper for P-40U, P-50U, P-60U, P-61U, P-63UM, P-65U, P-67UA, P-68U, P-90U, P-90W, P-91 series and CP-700DU/DSU printers. 215 prints/roll	610-064	4 rolls/box
KP-65HM b/w high density paper with blue background for P-40U, P-63UM, P-65U, P-67UA, P-68U, P-90U, P-90W, P-91 series, P-93W, P-93DW, CP-700UM, CP-700DU/DSU and CP-770DW printers. 215 prints/roll	610-314	4 rolls/box
KP-70S b/w super grade paper for P-70U, P-71U, P-75U, P-78U and P-500W printers, 115 prints/roll	610-065	4 rolls/box
KP-91HG high gloss thermal paper for P-91W, P-91DW, P-91DW/UB, P-93W and P-93DW printers.	610-812	4 rolls/box
CK-700 color paper for use with CP-700DU, CP-700U, CP-700DSU, CP-750U, CP-700UM, CP-710U and CP-770SW printers. For use with PK-700L and PK-700S ink ribbons. 200 prints/roll.	610-628	10 rolls/box
PK-700L large format printer ink ribbon for use with CP-700DU, CP-700U, CP-700DSU, CP-750U, CP-700UM, CP-710U and CP-770DW printers. 130 prints.	610-623	10 ribbons/box
PK-700S standard format printer ink ribbon for use with CP-700DU/DSU, CP-700U, CP-750U, CP-700UM, CP-710U, and CP-770DW printers. 200 prints.	610-624	10 ribbons/box
CK-800L color paper for CP-800UM & CP-800E series color printers. 160 prints, ink sheet included.	610-758	1 roll/box
CK-800S color paper for CP-800UM & CP-800E series color printers. 160 prints, ink sheet included.	610-759	1 roll/box
CK-900L standard paper, L-size ink sheet for CP-900UM, CP-900DW and CP-910U printers.	610-813	130 sheets/box
CK-900S standard paper, S-size ink sheet for CP-900UM, CP-900DW and CP-910U printers.	610-814	200 sheets/box
<b>MITSUBISHI PRINTERS</b>	<b>REF. #</b>	<b>QTY.</b>
NEW • Mitsubishi CP-30W color printer	610-1042	1
NEW • Mitsubishi CP-910U color thermal printer	610-1044	1

SONY PRINT MEDIA	REF. #	QTY.
UPP-110HA b/w dynamic range paper for UP-850, UP-880, UP-890MD, UP-D890 printers. 110mm x 18m. 240 prints/roll.	610-389	10 rolls/box
UPP-110HD b/w high density paper for UP-850, UP-D860, UP-870MD, UP-880, UP-890MD, UP-895MD, UP-D895, UP-897MD, UP-D897 printers. 110mm x 20m. 270 prints/roll.	610-072	10 rolls/box
UPP-110HG b/w high density gloss paper for UP-895MD, UP-D895, UP-897MD, UP-D897 printers. 110mm x 18m. 240 prints/roll.	610-754	10 rolls/box
UPP-110S/5 (replaces UPP-110S/4) b/w standard paper for UP-810, UP-811, UP-850, UP-860, UP-870, UP-880, UP-890MD, UP-895MD, UP-D895MD, UP-897MD, UP-D897 printers. 110mm x 20m. 210 prints/roll.	610-071	10 rolls/box
UPP-210HD b/w high density paper for UP-910, UP-930, UP-960, UP-980 printers. 210mm x 25m. 130 prints/roll.	610-619	5 rolls/box
UPP-210SE b/w standard paper for UP-910, UP-930, UP-960, UP-980 printers. 210mm x 25m. 130 prints/roll.	610-620	5 rolls/box
UPT-210BL blue transparency film for use with UP-980 printer. 210mm x 12.5m. 42 prints/roll.	610-621	5 rolls/box
UPC-21S color print pack (S-size paper) for the UP-20, UP-21MD, UP-D21MD, UP-D23MD printers. 3 ink ribbons and 1 cleaning sheet.	610-858	240 prints/box
UPC-21L color print pack (L-size paper) for the UP-20, UP-21MD, UP-D21, UP-D21MD, UP-D23MD printers. 4 ink ribbons, 1 cleaning sheet.	610-857	200 prints/box
UPC-55 color print pack for the UP-55MD, UP-D55MD printer. Includes 2 ink ribbons.	610-931	200 sheets/box
UPC-510 color print pack for UP-D50, UP-51MD, UP-51MDU printers. Includes ink ribbon.	610-757	200 sheets/box
UPC-1010/1 color print pack for UP-1200, UP-1800MD, UP-1850MD, UP-D1500CN printers. Ink cartridge included.	610-395	100 prints/box
UPC-2010/3 color print pack for the UP-2100, UP-2300, UP-D2500, UP-2900MD printers. Ink cartridge included.	610-617	200 prints/box
UPC-3010 color printing pack for UP-3000, UP-3030 printers. Ink cartridge included.	610-247	80 prints/box
UPC-5010A/2 color printing pack for UP-5000, UP-5000W, UP-5200MD printers. Ink cartridge included.	610-555	100 prints/box
UPC-5020A b/w printing pack for UP-5000, UP-5200MD printers. Ink cartridge included.	610-753	100 prints/box
UPC-5510 color printing pack for UP-5500, UP-D5500, UP-5600MD, UP-5600MDU printers. Ink ribbon included for re-loading into existing cartridges.	610-390	200 prints/box
SONY PRINTERS	REF. #	QTY.
SONY UP-D897MD medical grade digital black and white printer	610-995	1
SONY UP-897MD black & white printer	610-953	1
Digital Recorders	REF. #	QTY.
SONY DVO-1000MD medical grade DVD recorder	610-994	1

**FOR MORE  
INFORMATION...**

Search CIVCO's easy-to-use cross-reference charts to find compatible needle guides, transducer covers, ultrasound supplies and other accessories recommended for use with your Siemens Ultrasound System at:

**WWW.CIVCO.COM/SIEMENS**

*Needle Guidance*



*Transducer Positioning*



*In-Plane Needle Guidance*



*Infection Prevention*



For international distribution, visit CIVCO Worldwide at [WWW.CIVCO.COM](http://WWW.CIVCO.COM)

Worldwide Headquarters: 102 First Street South, Kalona, Iowa 52247-9589 USA  
319/656-4447, 800/445-6741 • Fax: 319/656-4451, 877/329-2482 • [WWW.CIVCO.COM](http://WWW.CIVCO.COM)

COPYRIGHT © 2011. ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. BRACHY BALLOON, CIV-FLEX, MULTI-PRO 2000, NEOGUARD, PROCOVERS, QUICK-MOUNT, SURGL, SURGI-BOOT, SURGI-TIP, ULTRA-PRO II AND ULTRA-PRO 3. ARE ALL TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE THE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. CIVCO PRODUCT NAMES MAY BE COVERED BY ONE OF THE FOLLOWING U.S. PATENT NOS. 5,941,889, 6,216,607 AND 7,909,815. PRINTED IN USA 2011R-2545